

KALOS

Clearing tracheobronchial secretions:

In order to remove the secretions in a non-invasive way, avoiding the use of suction catheters and reducing infection risk and discomfort, we developed a medical device with various techniques.

Cough assist:

The choice between manual, automatic and synchronized mode makes the therapy customizable. The possibility of save three therapy profile helps the patient in the daily treatment. Four insufflation air flow rates makes the therapy more comfortable.

Free aspire:

During expiration, the expired air is accelerated proportionally to the flow. This technique has been used for years in order to remove the secretions in a non-invasive way and without using positive or negative air pressure. Completely safe for the patients, it can be used as well with infant.

Continuous pressure:

Used if a continuous positive pressure is needed in order to make a lung expansion exercises.

Peak Flow Test

A fun interface allows to evaluate the estimated peak flow into three slots: under 160 Vmin, between 160 Vmin and 270 Vmin, more than 270 Vmin.





TECHNICAL DATA

Mode / Functions Free Aspire

Level adjustable from 1 to 5

Therapy time adjustable from 0 to 30min

(E₁₉₃₆

Cough Assist

Modes manual/auto /sync Inspiratory pressure from 5 to 70 cmH2O Expiratory pressure from -5 to -70 cmH2O Inspiratory time from 0.5 to 5 sec. Expiratory time from 0.5 to 5 sec. Pause time from 0.5 to 10 sec. Rise time adjustable from 1 to 4 from 5 to 20 Hz Pulsate frequency Pulsate amplitude from 1 to 5 cmH2O Duration of therapy adjustable

Continuous pressure

Pressure from 4 to 20 cmH2O

Duration of therapy adjustable from 0 to 15min

AC Voltage Source 100-240V AC - 50-60Hz
DC internal power source NiMH 3.2Ah 12V cc
Electrical Absorption 160VA (70W)
Type of Protection Against Electric Shock class II
Degree of protection against electric shock type BF

Degree of protection against electric shock type B

Degree of Protection against ingress of water IP21

Dimensions and weight 27 x 30 x 18 cm - 2,9 Kg

International standards:

EN 60601-1, EN 60601-1-11, EN 60601-1-8, EN 10651-6

Safety of Electromedical Instruments, General Requirements Electromagnetic Compatibility in accordance with EMC Directive. Medical Products Research srl declares that Kalos complies with the following EMC standard: EN60601-1-2 Illed.

This technical specifications is updated at the moment of printing and they can change depending on our politics of continuous improvement. ©MPR SrI



